

VITARI™

Scale without the premium.



MEET VITARI

High-throughput sequencing has been built on compromise. Scale or flexibility. Performance or usability. Efficiency or quality.

VITARI changes that. Scale and flexibility coexist. Performance is intuitive. Quality is uncompromised. It's built for what comes next.

- Scale output with **up to 10 billion reads or 3 TB** data per dual-sided run
- Lower costs at scale with a **\$100 genome**
- Rely on high-confidence data **>90% Q30 quality**
- Flexible, right-sized run configurations with **2 independent** flow cells with **6 lanes each**



~36 hour run times (2 x 150)
aligned with real lab workflows

High throughput
performance in a compact
benchtop platform

Simplify logistics with
ambient shipping of the
sequencing reagent cartridge

Compatible with Cloudbreak
Freestyle™ and Trinity™
2x150 and 2x75 kits



Shipping in H2 2026. Pre-order yours today:

Design with intent

At Element, we innovate to meet researchers where they are today and where their science is headed tomorrow. We build platforms around real lab needs: reliable performance, practical workflows, and room to grow so you can focus on discovery. And we keep evolving what we deliver, so your platform stays aligned with new applications, new questions, and new possibilities.

Explore our platforms:



VITARI

A high-throughput benchtop sequencing platform designed to remove compromises and is built for what comes next.



AVITI24™

Seamlessly integrates high-quality sequencing and onboard 5D multiomics into a single dual platform.



AVITI™

A flexible mid-throughput sequencing platform designed for everyday reliability and high-quality results across a range of applications.

Information in this document is provided for research use only and is subject to change without notice.

© 2026 Element Biosciences, Inc. All rights reserved. ABC, Adept, Avidite, Avidite Base Chemistry, AVITI, AVITI24, AVITI Dx, Catalyst, CB Freestyle, Cloudbreak, Cloudbreak Freestyle, CytoCanvas, Elembio, ElemBio, Elembio Catalyst, Elembio Cloud Custom Designer, ElemBio Cloud Custom Designer, Element Biosciences, Elevate, Teton, Teton Atlas, Trinity, Trinity Freestyle, UltraQ, and VITARI are trademarks and/or registered trademarks of Element Biosciences, Inc. in the United States and/or other countries. All other trademarks, service marks, and logos are the property of their respective owners. Use of third-party marks is for identification purposes only and does not imply endorsement, sponsorship, or affiliation.

Visit elementbiosciences.com for more information.