




Deep Well Plates are used for sample collection, in-vitro growth chambers and longterm storage. Made of virgin polypropylene, they offer excellent chemical resistance and will withstand temperatures down to -80°C.

- 1.1ml and 2.0ml capacities, in durable medical grade polypropylene
- 1.1ml features an ultra-low profile (27mm) for reduced space requirements
- Excellent chemical resistance
- Standard SLAS (SBS) footprint
- Centrifugation to 4,000 RCF
- 'V' bottom allows complete sample recovery
- Ultra-flat surface ensures proper sealing of heat film



DWP-11-RDU


CATALOG NO.	VOLUME TO BRIM	WORKING VOLUME	VOLUME WITH MAT	DESCRIPTION	PACK	CASE
DWP-20-SQV	2.3ml	2.2ml	1.9ml	96-well Deep Well Plate, Polypropylene, 2.0ml, Clear, V-Bottom, Square Wells	5 Plates	10 Packs
DWP-20-SQU	2.2ml	2.1ml	1.8ml	96-well Deep Well Plate, Polypropylene, 2.0ml, Clear, Round Bottom, Square Wells	5 Plates	10 Packs
DWP-20-RDU	1.9ml	1.8ml	1.5ml	96-well Deep Well Plate, Polypropylene, 2.0ml, Clear, Round Bottom, Round Well	5 Plates	10 Packs
DWP-11-RDU	1.2ml	1.1ml	1.0ml	96-well Deep Well Plate, Polypropylene, 1.1ml, Clear, Round Bottom, Round Wells	5 Plates	10 Packs




Silicone Sealing mats are designed for **high throughput screening** and **storage applications**. Made of research grade silicone, they can be autoclaved and following a bleach wash and ethanol rinse protocol, they can also be reused. Silicone Sealing mats, pierceable and self-sealing, are ideal for automated processes and can also be used in subzero environments.




DWP-20-SQV



DWP-20-RDU



Round Bottom V-Bottom



Aluminum Film

Polyester Film

Sealing mats & films

CATALOG NO.	DESCRIPTION
SM-2ML-RD	Sealing mat, for 1.1 and 2ml 96-well plates, Round wells, Research Grade Silicone, Pierceable
SM-2ML-SQ	Sealing mat, for 2ml 96-well plates, Square wells, Research Grade Silicone, Pierceable
SM-750UL-RD	Sealing mat, for Greiner 96-well Masterblock plates, Round wells, Research Grade Silicone, Pierceable
SF-AS-400	Aluminum sealing film for PCR and Cold Storage
SF-TS-600	Polyester sealing film for PCR and Storage
SF-TS-200	General application sealing film
SM96-SEPTA3100	SeptaMat for ABI Multi-Capillary Sequencing Instruments

Packaging Description:

Mats - 10 mats/pack, 5 packs/case
Films - 100 films/pack, 5 packs/case

