

UNIVERSAL FIT TIPS

Molded in virgin medical grade polypropylene, NOVASBIO Pipette Tips provide exceptional performance especially when pipetting microvolumes.

Universal fit - compatible with most popular single and multichannel pipettes

Available in racked and refill formats

The thin-wall design creates a perfect seal and excellent compatibility with most popular brands of pipettors. These universal-fitting pipette tips are engineered to fit practically all single and multichannel pipettes.

NOVASBIO pipette tips can be autoclaved.



EXTENDED LENGTH MICROTIPS REACH DOWN TO THE BOTTOM OF TUBES, REDUCING THE RISK OF CONTAMINATION

REFILL SYSTEM

REUSE, REDUCE AND RECYCLE!

Available in three popular volumes, 10µl, 200µl and 1000µl, NOVASBIO Refill System utilizes minimum packaging and so reduces waste and environmental concerns

Uses 42-55% less plastics in the packaging compared to standard racked tips

All materials used for racks and packaging are recyclable

Kraft paper boxes, 100% FSC certified to ensure the packaging boxes come from responsibly managed forests

Minimized printing/labels which are also recyclable



CATALOG NO.	VOLUME	COLOR	DESCRIPTION	TIP LENGTH	PACK	CASE
T-10-C-R	10µl	Clear	Racked, Graduations, Non-Sterile	31mm	10 x Rack/96 Tips	5 Packs
T-10-C-R-S	10µl*	Clear	Racked, Graduations, Pre-Sterilized	31mm	10 x Rack/96 Tips	5 Packs
T-10-C-L-R-S	10µl*	Clear	Low Retention, Racked, Graduations, Pre-Sterilized	31mm	10 x Rack/96 Tips	5 Packs
RFL-T10	10µl	Clear	Refill Racked, Graduations, Non-Sterile	31mm	12 x Rack/96 Tips	5 Packs
T-10XL-C-R	20µl	Clear	Extended, Racked, Non-Sterile	45mm	10 x Rack/96 Tips	5 Packs
T-200-C-R	200µl	Clear	Racked, Non-Sterile	52mm	10 x Rack/96 Tips	5 Packs
T-200-C-R-S	200µl	Clear	Racked, Pre-Sterilized	52mm	10 x Rack/96 Tips	5 Packs
T-200-C-L-R-S	200µl	Clear	Low Retention, Racked, Pre-Sterilized	52mm	10 x Rack/96 Tips	5 Packs
RFL-T200	200µl	Clear	Refill Racked, Non-Sterile	52mm	10 x Rack/96 Tips	5 Packs
T-1000-C-R	1000µl	Clear	Racked, Non-Sterile	78mm	10 x Rack/100 Tips	5 Packs
T-1000-C-R-S	1000µl	Clear	Racked, Pre-Sterilized	78mm	10 x Rack/100 Tips	5 Packs
T-1000-C-L-R-S	1000µl	Clear	Low Retention, Racked, Pre-Sterilized	78mm	10 x Rack/100 Tips	5 Packs
T-1000XL-C-R	1000µl*	Clear	Extended, Racked, Graduations, Non-Sterile	83mm	10 x Rack/96 Tips	5 Packs
T-1000XL-C-L-R-S	1000µl*	Clear	Extended, Low Retention, Racked, Graduations, Pre-Sterilized	83mm	10 x Rack/96 Tips	5 Packs
RFL-T1000XL	1000µl*	Clear	Extended, Refill Racked, Graduations, Non-Sterile	83mm	8 x Rack/96 Tips	5 Packs

Available in both Bulk and Sterile versions upon request

* Full working volume 1200µl

